



آغا خان یونیورسٹی ہسپتال

The Aga Khan University Hospital
Clinical Laboratories

Specimen Collection and Transportation Information

Test Name	Mnemonic	Section	Specimen Type
Stool Concentration Detail Report	SDR	Clinical Microbiology	Fresh Stool

- **Preparation of the Patient:** NA
- **Type of Collection Container:** Wide mouth D/R container
- **Types and Amounts of Preservatives or Anticoagulants:** sample in formalin and PVA also accepted
- **Quantity of Specimen to be Collected:** 5ml of diarrheal stool, 1g of material or a walnut-sized portion of stool
 - **Minimum Quantity:** As above
- **Transportation Temperature:** Room temperature
- **Need for Special Timing for Collection (where applicable):** NA
- **Specimen Stability Information:** If delay in between 2-4hrs at Room temp, If delay more than 4 hrs. at refrigerated temp
- **Special Instruction:** NA
- **Rejection Criteria:** NA
- **Need for Appropriate Clinical Data, When Indicated (Patient history):** Required

Methodology	Reporting Scheme	Day Performed	Cut Off Time Main Lab	Cut Off Time Karachi Points	Cut Off Time Out of Karachi Points
Microscopy (Manual)	Within 24 Hours	Daily	00.01	10.00	Contact nearest Collection Point