



آغا خان یونیورسٹی ہسپتال

The Aga Khan University Hospital

Clinical Laboratories

Specimen Collection and Transportation Information

Test Name	Mnemonic	Section	Specimen Type
MEASELS IgM	MAB	Clinical Microbiology	Serum or Clotted Blood

- **Preparation of the Patient:** NA
- **Type of Collection Container:** Gel tube
- **Types and Amounts of Preservatives or Anticoagulants:** NA
- **Quantity of Specimen to be Collected:** 3-5cc
 - **Minimum Quantity:** 3cc
- **Transportation Temperature:** Room temp
- **Need for Special Timing for Collection (where applicable):** NA
- **Specimen Stability Information:** If delay in between 2-4hrs at Room temp, If delay more than 2 hrs at refrigerated temp
- **Special Instruction:** NA
- **Rejection Criteria:** Hemolyzed blood not accepted
- **Need for Appropriate Clinical Data, When Indicated (Patient history):** Not Required

Methodology	Reporting Scheme	Day Performed	Cut Off Time Main Lab	Cut Off Time Karachi Points	Cut Off Time Out of Karachi Points
ELISA	Done on 1st & 3rd Thu, Reported Next Day	Weekly	00.01	10.00	Contact nearest Collection Point